

CERNER TRAINING MANUALS

ACCESSION RESULT ENTRY



CONTENTS

LIST OF TABLES	v		
INTRODUCTION	1		
About this Procedure	1		
Accession Result Entry Modes	1		
CUSTOMIZING ACCESSION RESULT ENTRY	3		
Setup Result View	3		
Assay Sorting	4		
CUSTOMIZING DEMOGRAPHICS	5		
Helpful Tips	5		
Changing Demographics	5		
Restore to Defaults	8		
ENTERING RESULTS	9		
Manual Result Entry	9		
Result Types	11		
Entering Comments	11		
INTERPRETING RESULTS	13		
Flags	13		
Review Range	13		
Critical Results	13		
Linearities	14		
Reference Ranges	15		
VERIFYING RESULTS	17		
Verifying Critical Results	17		
DIFFERENTIAL MODE	19		
Orders with Differentials	19		
		Body Fluid Differentials	19
		CBC Differentials	24
		CORRECTION MODE	29
		Correcting Results	30
		INSTRUMENT QUEUE MODE	33
		Opening Instrument Queue	33
		Cosmetic Differences	34
		Using Instrument Queue	34
		ADDITIONAL INFORMATION	37
		Getting Assay Information	37
		Previous Results	38
		TROUBLESHOOTING	41
		Too Many Digits	41
		Decimals Not Allowed	41
		Unable to Enter Results	42
		ROUTING	45
		What is Routing?	45
		How Does It Work?	45
		Brackenridge	46
		Dell Children's	47
		Edgar B. Davis	48
		Hays	49
		Highland Lakes	50
		Northwest	51
		Southwest	52
		Williamson	53
		Seton Medical Center	54

LIST OF TABLES

1	Result Types	11	3	Manual Differentials and Codes	19
2	Result Flags	13	4	CBC Orders	24
			5	CBC Reflexing Rules	25

INTRODUCTION

Accession Result Entry (ARE) is the Cerner application used to enter results. It has several modes which allow for manual result entry, resulting from instruments, or performing differentials.

ARE can also show previous results, reference ranges, linearities, and critical ranges.

ABOUT THIS PROCEDURE

Since resulting is the same in most of the ARE modes,¹ CHAPTER INTERPRETING RESULTS pg.13 will cover resulting orders using Accession Result Entry.

¹ Differential Mode is the one exception.

The chapters for each mode will simply walk you through the differences.

ACCESSION RESULT ENTRY MODES

Depending on which tests are going to be resulted, ARE has different modes to make resulting easier.

ACCESSION MODE: can be used to enter results manually, or verify results performed by instruments.

DIFFERENTIAL MODE: is used to perform DIFFERENTIALS using Cerner. Mostly, it will be used for BODY FLUIDS.

It can be used during WAM DOWNTIMES to perform CBC differentials.

INSTRUMENT MODE: is used to review results which cannot be AUTO-VERIFIED by the instruments.



CORRECTION MODE: is used to correct results which have already been VERIFIED.

CUSTOMIZING ACCESSION RESULT ENTRY


SETUP RESULT VIEW

The RESULTS SPREADSHEET of ARE can be modified to display additional information. Specifically, the REFERENCE RANGE column² should be added.

² It's really handy.


 *Open* ARE by clicking the  icon from the App-Bar.

 *Click* **VIEW** from the menu bar.



 *Click* **CUSTOMIZE...**


 *Click* **RESULT DISPLAY...**



 *Check* ☒ **REFERENCE RANGE.**

 *Check* ☒ **SERVICE RESOURCE DISPLAY.**³

³ This option will add a column to ARE. When instruments perform results, the Instrument ID will be present in this column.



 **Un-Check** ☐ the other options.⁴

⁴ Or leave them Un-checked. They really don't add much.


Now, NEXT TO THE sample status is a column which lists the reference ranges for each assay.


Result:		Hierarchy:					
Procedure		Result	Flags	Status	Reference Range	Service ...	Time
Glucose Level	<input checked="" type="checkbox"/>			Pending	70 - 110		
BUN	<input type="checkbox"/>			Pending	6 - 20		
Creatinine	<input type="checkbox"/>			Pending	0.5 - 1.2	15.0 H	61 days
Sodium Level	<input type="checkbox"/>			Pending	136 - 145	136	47 days
Potassium Level	<input type="checkbox"/>			Pending	3.6 - 5.0	5.0	47 days
Chloride	<input type="checkbox"/>			Pending	101 - 111	101	47 days
CO2	<input type="checkbox"/>			Pending	21 - 31	28	47 days
AGAP	<input type="checkbox"/>			Pending	4 - 12	7	47 days
Calcium	<input type="checkbox"/>			Pending	8.5 - 10.5		

ASSAY SORTING

By default, the assays are sorted by ORDERABLE. This sorting method is very inconvenient if the assays are performed in multiple departments.

Sorting by SERVICE RESOURCE will group assays based on where they're performed.


 **Click VIEW** from the menu bar.



 **Click CUSTOMIZE...**

 **Click SORT ASSAYS.**

 **Click BY SERVICE RESOURCE.**



Now, ASSAYS ARE WILL be sorted by BENCH, or testing location, making ARE much easier to use.

CUSTOMIZING DEMOGRAPHICS

There are a few applications within Cerner⁵ which display demographic information.

⁵ ARE and Department Order Entry (DOE) are the two most common applications with this feature.

The information displayed can be modified to display any necessary information. This can be very useful if you find yourself searching for specific patient information.

Demographics		
Name: TEST, VON T...	Admitting dx:	Loc: Medical Center / ICUA SMC / ICUA / 01 / Se...
FIN: 6038342000	MRN:	Ordered by: TestIS, Z VMiller
DOB: 12/10/1984	Gender: Male	Admitted: 5/15/2015 9:34
Order Comment:		
Location: ICUA SMC / ICUA / 01		

HELPFUL TIPS

Before you begin, here are a few tips you should keep in mind if you decide to create your own demographics setup.


FOLLOW THESE TIPS for the best results.

- **NUMBER OF COLUMNS:** The fewer the columns the better.⁶
- **AVOID TMI:** Only include fields you *know* you'll need.
- **SPLIT LOCATIONS:** Separate the patients floor from the room number.

⁶ Fewer columns means information will not be cutoff.


CHANGING DEMOGRAPHICS

 Click **VIEW** from the menu bar.




 Click **CUSTOMIZE...**

 Click **DEMOGRAPHICS...**




THIS WINDOW displays the VALUES currently set in each of the demographics FIELDS.




 Set the columns count to 2.⁷


⁷ It's located on the top left corner of the window.




CERNER WILL ASK if you're sure.



 Click **Yes**

NEXT TO EACH FIELD in the DEMOGRAPHICS DATA ELEMENTS box is an ellipses  button. Clicking it will open the SET DEMOGRAPHIC DATA ELEMENTS WINDOW.

The two important fields are VALUE⁸ and LABEL.⁹



⁸ VALUE is what information will be displayed.

⁹ LABEL is the name it is given. This can be customized.


FOR EACH OF THE 10 fields listed:



 Click the  button.

 Set the VALUES to match the image below.

 Modify the LABEL, if needed.¹⁰



¹⁰ Sometimes it's helpful to make it shorter.



 Click  when finished.

RESTORE TO DEFAULTS

At any the values can be restored to their original settings.

 Click the  button.

ENTERING RESULTS

This section will describe the process of **MANUAL RESULT ENTRY**, **PERFORMING RESULTS**, and **VERIFYING RESULTS**.

In addition, it describes **RESULT FLAGS**, **REFERENCE RANGES**, **CRITICAL RANGES**, **LINEARITIES**,

ACCESSION MODE is the default resulting mode for ARE. It can be used to enter results manually, or verify results performed by instruments.

MANUAL RESULT ENTRY

Open ARE by clicking the from the App-Bar.

Enter the accession number.¹¹

Click

¹¹ See the **ACCESSION NUMBERS** procedure for tips on entering Accession Numbers.

PathNet General Lab: Accession Result Entry

Task Mode View Help

Selection
 Accession: 15-292-000102 Retrieve
 Procedure: All RT - Routine
 Test Site: All


Demographics
 Na... TEST, V... Admitting dx: Loc: Medical Center / ICUA SMC / IC...
 FIN: 60383420... MRN: Ordered by: TestIS, Z VMiller
 DOB: 12/10/19... Gender: Male Admitted: 5/15/2015 9:34
 Order Comment:
 Location: ICUA SMC / ICUA / 01

Result: Hierarchy:


Procedure	Result	Flags	Status	Prev	Time
Glucose Level			Pending		
BUN			Pending		
Creatinine			Pending	15.0 H	61 days
Sodium Level			Pending	136	47 days
Potassium Level			Pending	5.0	47 days
Chloride			Pending	101	47 days
CO2			Pending	28	47 days
AGAP			Pending	7	47 days
Calcium			Pending		

Perform... Verify...

Order: BMPNL; Status: Ordered; Procedure: Glucose Level; Location: SMC Remisol; Collected: 10/19/2015 20:52 A3291 TPERR005 10/19/2015 21:15



 *Click* on the RESULT CELL to enter a result.

Procedure	Result
Glucose Level	
BUN	
Creatinine	


 *Enter* the result.¹²

¹² Cerner will fill any missing decimal points with ZEROS.

Procedure	Result
Glucose Level	98
BUN	
Creatinine	

 *Hit*  ENTER on the keyboard to move to the next result.

Procedure	Result
Glucose Level	98
BUN	
Creatinine	

 *Continue* entering results.

Procedure	Result	Flags	Status	Reference Range
Glucose Level	98		Pending	70 - 110
BUN	7		Pending	6 - 20
Creatinine	1.0		Pending	0.5 - 1.2
Sodium Level	136		Pending	136 - 145
Potassium Level	5.0		Pending	3.6 - 5.0
Chloride	101		Pending	101 - 111
CO2	28		Pending	21 - 31
AGAP	7		Pending	4 - 12
Calcium	10.0		Pending	8.5 - 10.5

RESULT TYPES

Every assay has a default RESULT TYPE. This is the result type which will be used most often to result the assay.

In situations where it's not possible to use the default,¹³ the RESULT TYPE will need to be converted to an appropriate alternative.

¹³ e.g. Unable to perform calculations.

RESULT TYPES


There are several ways results can be entered using ARE:

TYPE	MEANING	HOW IT'S ENTERED
ALPHA	Pre-composed results.	Select via Drop down
CALCULATION	Calculated values.	Automatically Entered.
FREETEXT	Short text result.	Type result.
NUMERIC	Measured values.	Type result.
TEXT	Long text results.	Word Processor.


Table 1: Result Types

CHANGING RESULT TYPES


There are situations where it may be necessary to change a result type. For instance: If a TEG doesn't split, the result needs to be changed to an **ALPHA RESPONSE**.

 *Right Click* on the result cell.

 *Select* **CONVERT RESULT**.


 *Select* the appropriate RESULT TYPE.¹⁴


¹⁴ If a RESULT TYPE for that assay is unavailable it will be disabled in this menu.



ENTERING COMMENTS

Additional information can be attached to individual results or entire orders by adding Comments and/or notes.


info: For a detailed description on comments:  Refer to the **COMMENTS** documentation.

 *Select* an assay which needs a comment.¹⁵

 *Click*  from the TOOLBAR



THE COMMENTS WINDOW will open.



ORDER COMMENT: A comment applied to the entire order.¹⁶

ORDER NOTE: An internal note applied to the entire order.

RESULT COMMENT: A comment attached to a specific result.¹⁷

RESULT NOTE: An internal note applied to a specific result.

¹⁵ If the comment or note will be attached to the entire order, then simply select any result.

¹⁶ In this example it is applied to the whole Basic Metabolic Panel.

¹⁷ In this example it is applied to just the Creatinine result.

INTERPRETING RESULTS

FLAGS

Any result which *may* require special attention will be flagged. These flags appear as letters in the FLAGS Column of ARE.

TYPE	SYMB.	MEANING
REVIEW	R	The results needs to be reviewed before verifying.
DELTA	D	The result is has changed significantly.
HIGH	H	The result is above the Normal Range.
LOW	L	The result is below the Normal Range.
CRITICAL	C	The result is Critical.
LINEARITY	N	The result is outside of the defined Linearity.
NOTIFY	@	The result needs to be called to the floor.
INTERPRETATION	(t)	The result is an an interpretation.
COMMENT	f	The result contains a Comment or Note.

Table 2: Result Flags

THE EXAMPLE BELOW shows multiple flags¹⁸ for the result.

Procedure		Result	Flags	Status	Reference Range
Glucose Level	<input checked="" type="checkbox"/>	> 2400	HCRNf	Pending	70 - 110

¹⁸ Specifically:

- HIGH
- CRITICAL
- REVIEW
- LINEARITY
- COMMENT

REVIEW RANGE

The REVIEW RANGE is used to prevent results from auto-verifying, and to alert the user that action might need to be taken.

CRITICAL RESULTS

ARE will emphasize critical results by changing the color to RED and adding the CRITICAL FLAG (C).

Sodium Level	<input checked="" type="checkbox"/>	190	HCR	Pending	136
Potassium Level	<input checked="" type="checkbox"/>	2.0	LCR	Pending	5.0

LINEARITIES

The LINEARITIES have been defined as the minimum and maximum reportable values **AFTER** dilution.

Any result that does not fall within the LINEAR RANGE will automatically be converted to its inequality,¹⁹ and the LINEARITY FLAG (N)²⁰ will be added.

info: The best way to manually enter results which are out of linearity is to enter a result which is outside of the linearity.²¹

When the instrument sends over results which are outside of the LINEAR RANGE, it will automatically be converted to its inequality.

LINEARITIES AND CALCULATIONS

When a component of a calculation falls outside of its LINEAR RANGE, Cerner will perform the calculation using the linear limit of the assay.²²

Finally, the appropriate inequality will be attached to the calculated result.

THE FOLLOWING TWO EXAMPLES show what happens if a component is above, and below linearity.

BELOW LINEARITY:

Procedure		Result	Flags
Hours Collected	<input checked="" type="checkbox"/>	24	
Total Volume	<input checked="" type="checkbox"/>	3265	
U Timed Protein	<input checked="" type="checkbox"/>	< 6	N
U24 Protein	<input checked="" type="checkbox"/>	< 196	H

¹⁹ e.g. A glucose of 2,401mg/dL will convert to >2,400mg/dL.

²⁰ N is used because L was taken by the Low flag.

²¹ Protip brought to you by Cpt. Obvious.

²² e.g. If the U TIMED PROTEIN IS <6MG/DL, CERNER WILL PERFORM THE CALCULATION USING A VALUE OF 6MG/DL.


In this case, we can accurately say the U24 Protein is <196mg/day.

ABOVE LINEARITY:

Procedure		Result	Flags	Status
Hours Collected	<input checked="" type="checkbox"/>	24		Pending
Total Volume	<input checked="" type="checkbox"/>	3265		Pending
U Timed Protein	<input checked="" type="checkbox"/>	> 1500	N	Pending
U24 Protein	<input checked="" type="checkbox"/>	> 48975	H	Pending


In this case, we can accurately say the U24 Protein is >48,975mg/day.

IF MULTIPLE COMPONENTS of the calculation are outside of LINEARITY, Cerner will not be able to perform the calculation. Instead, a message will appear stating it's unable to perform the calculation.



IN THIS SITUATION the calculated result needs to be manually entered as UNABLE TO CALCULATE.²³

²³ For more information [Refer to: RESULT TYPES pg. 11.](#)

 **Convert** The result to an ALPHA RESPONSE.

 **Select** UNABLE TO CALCULATE.

info: If you're unable to convert it to an ALPHA RESPONSE, convert it to a FREETEXT and manually type: *"Unable to Calculate."*

LINEARITIES AND DELTA CHECKS

Cerner will **NOT** be able to perform DELTA CHECKING if any of the two values are outside of linearity.

REFERENCE RANGES

Results which are lower than the defined REFERENCE RANGE will be flagged as Low (L), results which are higher will flag as HIGH (H.)²⁴

²⁴ What else can I say?

VERIFYING RESULTS

IF YOU'D LIKE TO SAVE the results, but **NOT** release them:

 Click  to PERFORM²⁵ the results.

²⁵ This will save the results; they can be modified or VERIFIED at a later time.



IF THE RESULTS ARE GOOD to be released:


 Click  to FINALIZE²⁶ the results.

²⁶ This will post the results to the patient's chart, they cannot be modified without CORRECTIONS.


VERIFYING CRITICAL RESULTS

ARE will do it's best to make sure that a comment has been placed on the result before allowing the result to be verified.


WHEN THE  or  button is pressed **AND** there is no COMMENT FLAG ARE will display one of the following alerts.




This is displayed if there is only one critical value.

CLICKING THE  button will open a COMMENTS WINDOW.²⁷

²⁷ **info:** The comment entered into this window will be applied to the CRITICAL RESULT.




This is displayed if there are multiple critical values.




CLICKING THE  button will open a BATCH COMMENTS WINDOW.²⁸

²⁸ **info:** The comment entered into this window will be applied to every CRITICAL RESULT.


⚠ WARNING: ARE will only check for the presence of a comment, it has no way of knowing if the comment is a **Critical Call** comment.



Quick Reminder:

- 🖨 Type **CALRED**
- 🖨 hit  **F9** to expand.
- 🖨 hit  **F3** to move.
- 🖨 hit  **F5** for date.

info: The great thing about this feature is that it's a quick way to apply the CRITICAL CALL COMMENT to every critical result.

For more information on adding comments  refer to the **COMMENTS** documentation.

DIFFERENTIAL MODE

Resulting differentials in Cerner is much different than resulting anything else. For that reason, Accession Result Entry has a mode specifically for performing them.

Differential mode will perform the normalization calculations and with CBC's it will calculate the absolute values.

⚠ IMPORTANT: All manual differentials should be entered using differential mode!!!²⁹

²⁹ If I could make this statement flash I would.

ORDERS WITH DIFFERENTIALS

CBCs: Manual differentials will be performed using WAM.³⁰

BODY FLUIDS: Differentials will be reflexed, as needed, after the body fluid counts have been VERIFIED.

³⁰ If WAM is down: [Refer to the DOWNTIME PROCEDURE](#) for information on adding Manual Differentials.

HERE IS A LIST of possible differential types:

FLUID	DIFFERENTIAL NAME	CODE
BODY FLUID	Body Fluid Differential	FLDIF
BAL FLUID	BAL Differential	BALDIF
CBC	Manual Differential	MDIFF
CSF	CSF Differential	CSFDIF
JOINT FLUID	Joint Fluid Differential	JNTDIF
PERICARDIAL FLUID	Pericardial Differential	PCFDIF
PERITONEAL FLUID	Peritoneal Differential	PERDIF
PLEURAL FLUID	Pleural Fluid Differential	PLFDIF

Table 3: Manual Differentials and Codes

HOW CERNER DETERMINES if a MANUAL DIFFERENTIAL is needed:


BODY FLUIDS: The result of the **WBC COUNT** is >5/ μ L.

CBC: The result of the **DIFFERENTIAL** assay.



BODY FLUID DIFFERENTIALS

In order for a manual differential reflex on a body fluid, the COUNT needs to be entered and VERIFIED. This is done using ARE IN ACCESSION MODE.

REFLEXING MANUAL DIFFERENTIALS


 *Enter* the results for the body fluid count.

Result:				
Procedure		Result	Flags	Status
Source Pleural Fluid	<input checked="" type="checkbox"/>	Left Pleur. Fl		Performed
Color Pleural Fluid	<input checked="" type="checkbox"/>	Straw		Pending
Character Pleural Fluid	<input checked="" type="checkbox"/>	Cloudy		Pending
WBC Pleural Fluid	<input checked="" type="checkbox"/>	12		Pending
RBC Pleural Fluid	<input checked="" type="checkbox"/>	230		Pending

 *Click*  when finished.³¹

³¹ The Differential will not reflex until the results are VERIFIED.

IF A MANUAL DIFFERENTIAL is added, Cerner will present a **DISCERN NOTIFICATION**.



CHANGING TO DIFFERENTIAL MODE

 *Open* the ACCESSION NUMBER in ARE.³²


³² The differential components should now appear in the ARE window.

Result:		Hierarchy:					
Procedure		Result	Flags	Status	Reference Range	Service Resource Disp...	Prev Time
Source Pleural Fluid	<input checked="" type="checkbox"/>	Left Pleur. Fl		Verified		SMC Heme Man	Right Pleural Fluid 5 months
Color Pleural Fluid	<input checked="" type="checkbox"/>	Straw		Verified		SMC Heme Man	Colorless 5 months
Character Pleural Fluid	<input checked="" type="checkbox"/>	Cloudy		Verified		SMC Heme Man	Hazy 5 months
WBC Pleural Fluid	<input checked="" type="checkbox"/>	12		Verified		SMC Heme Man	50 5 months
RBC Pleural Fluid	<input checked="" type="checkbox"/>	230		Verified		SMC Heme Man	50 5 months
Segs Pleural Fluid	<input checked="" type="checkbox"/>						6 13 days
Band Pleural Fluid	<input checked="" type="checkbox"/>						2 13 days
Lymph Pleural Fluid	<input checked="" type="checkbox"/>						80 13 days
Mono Histo Pleural Fluid	<input checked="" type="checkbox"/>						4 13 days
Eos Pleural Fluid	<input checked="" type="checkbox"/>						6 13 days
Basophil Pleural Fluid	<input checked="" type="checkbox"/>						2 13 days
Meso Pleural Fluid	<input checked="" type="checkbox"/>						
Other Cells Pleural Fluid	<input checked="" type="checkbox"/>						
Total Cells Counted Pleural Fluid	<input checked="" type="checkbox"/>	0		Pending			
Comment Pleural Fluid	<input checked="" type="checkbox"/>						


 **IMPORTANT:** Do NOT enter the results here!

 *Click* **MODE** from the menu bar.

Click **DIFFERENTIAL...**




THE SELECT ACCESSION box will display:



Click the **...** button.

Select the DIFFERENTIAL PROCEDURE from the list.³³


³³ If a differential is not there, it means the count has not been VERIFIED



Click **OK**

Select the **OPTION** using the dropdown menu.³⁴


³⁴ There should only be one option here.



Click **OK**


THE DIFFERENTIAL mode will open.

ABOUT DIFFERENTIAL MODE




SELECTION BOX: Allows you to change the number of counting Cells.³⁵

³⁵ For body fluids this will default to 300 Cells.



DIFFERENTIAL BOX: This displays the cells counted.




CELLS COUNTED: Displays the number of Cells counted. When finished, it will normalize the count to a percentage.


HOTKEY: The Key used to add the cell type.









PROCEDURE: A list of available cell types for the differential.

PERFORMING THE DIFFERENTIAL

The DIFFERENTIALS are performed using the 10-key pad.³⁶

³⁶ The  keys for each cell can be found immediately to the right of its name e. g. Segs are 1.

 *Count* the cells until you're finished.


Procedure		Result
Segs Pleural Fluid	1 	196
Lymph Pleural Fluid	2 	14
Mono Histios Pleural Fluid	3 	90
Band Pleural Fluid	4 	
Eos Pleural Fluid	5 	
Basophil Pleural Fluid	6 	
Meso Pleural Fluid	7 	
Other Cells Pleural Fluid	9 	


info: If you've run out of cells to count:

Click the STOP COUNTING icon  on the TOOL BAR

WHEN THE TOTAL CELLS counted reaches the set value,³⁷ an alert will appear.









³⁷ In this case, it's 300.





 *Click*  to acknowledge.³⁸

³⁸ This will normalize the differential.

NOTICE THAT THE CELLS have been normalized.

Procedure		Result
Segs Pleural Fluid	1 	65
Lymph Pleural Fluid	2 	5
Mono Histios Pleural Fluid	3 	30
Band Pleural Fluid	4 	
Eos Pleural Fluid	5 	
Basophil Pleural Fluid	6 	
Meso Pleural Fluid	7 	
Other Cells Pleural Fluid	9 	

OPTIONAL:

 *Click*  to save the results.

FINISHING UP


Before the results can be verified, the total number of cells counted needs to be entered.


info: Unfortunately, Cerner is unable to keep track of this information.

To DO THIS, ARE needs to be in **ACCESSION MODE**



 Click **MODE** from the menu bar.

 Click **ACCESSION**.



 Enter the NUMBER OF CELLS COUNTED in the TOTAL CELLS COUNTED result cell.


Total Cells Counted Pleural Fluid	<input checked="" type="checkbox"/>	300		Pending
Comment Pleural Fluid	<input type="checkbox"/>			

 Click  to verify the results

CBC DIFFERENTIALS

! IMPORTANT: CBC Differentials should be performed in WAM. The one exception is during **WAM DOWNTIMES**.

How CBCs Work

In Cerner, the components of a CBC are treated as separate orders. Cerner and WAM will do all the  magic required to switch between **MANUAL** and **AUTOMATED DIFFERENTIALS**.³⁹

³⁹ **info:** The ADIFF won't be added until the CBC is LOGGED-IN to the laboratory.

ORDER	DESCRIPTION
HEMGRM	Hemogram
CBC	Hemogram portion of CBC
ADIFF	Automated Differential
MORPHOLOGY	Slide Morphology
MDIFF	Manual Differential and Morphology

Table 4: CBC Orders

MANUALLY REFLEXING

⚠ IMPORTANT: Again, this only needs to be done if WAM is down.

CBCs have a result called **DIFFERENTIAL?**. This is used to trigger the reflexing of MORPHOLOGY or a MANUAL DIFFERENTIAL WITH MORPHOLOGY.

This table describes what happens with each of the **DIFFERENTIAL?** options.

SELECTION	CANCELED	ADDED	FINAL ORDER
(BLANK)	-	-	CBC, ADIFF
AUTO+MORPH	-	MORPHOLOGY	CBC, ADIFF, Morphology
MANUAL	ADIFF	MDIFF	CBC, MDIFF

Table 5: CBC Reflexing Rules

Select the appropriate option.⁴⁰

⁴⁰ Leave Blank if it's an AUTO DIFFERENTIAL without morphology.

Click

THIS WILL TRIGGER the reflexing rules. However, we won't see the changes until ARE is refreshed.

Click the ACCESSION FIELD.⁴¹

⁴¹ ARE will clear.



THE RESULTS SHOULD now reflect the changes.⁴²


If an AUTOMATED DIFFERENTIAL + MORPHOLOGY was selected, the MORPHOLOGY assays should appear.

⁴² If they don't, try again. Sometimes it takes a few seconds for the changes to be made.


PERFORMING MANUAL DIFFERENTIALS

Click **MODE** from the menu bar.

Click **DIFFERENTIAL...**




THE SELECT ACCESSION box will display:



Click the button.


Select **MDIFF**




Click OK

 *Select* the **OPTION** using the drop down menu.⁴³

⁴³ There should only be one option here.



 *Click* 



CBC MANUAL DIFFERENTIALS are performed similarly to Body fluid Differentials. The keys are set up to match WAM.

The right hand side has a list of morphology options; they can be entered using the DROP DOWN lists.


CORRECTION MODE


CORRECTION MODE is used to modify results, or result comments of assays which have already been VERIFIED.

When results are corrected their flags will be re-evaluated, and calculations⁴⁴ which use this value will be recalculated.


⁴⁴ Some calculations are more complicated and may need to be corrected manually. Examples are **PLATELET MAPPING, AND LIPIDS.**

 Click **MODE** from the menu bar.


 Click **CORRECTION.**




THE CORRECTION MODE ICON  will indicate that ARE is in CORRECTION MODE.



CORRECTING RESULTS

 **Select** the result which needs to be modified.


Procedure		Result	Flags	Status	Reference Range
Glucose Level	<input type="checkbox"/>	110		Verified	70 - 110
BUN	<input type="checkbox"/>	20		Verified	6 - 20
Creatinine	<input type="checkbox"/>	1.0		Verified	0.5 - 1.2
Sodium Level	<input checked="" type="checkbox"/>	150	H	Verified	136 - 145
Potassium Level	<input type="checkbox"/>	5.0		Verified	3.6 - 5.0
Chloride	<input type="checkbox"/>	101		Verified	101 - 111
CO2	<input type="checkbox"/>	29		Verified	21 - 31
AGAP	<input type="checkbox"/>	20	H	Verified	4 - 12
Calcium	<input type="checkbox"/>	8.5		Verified	8.5 - 10.5

 **Enter** the new result.


 **Hit**  **ENTER** on the keyboard.

Procedure		Result	Flags	Status	Reference Range
Glucose Level	<input type="checkbox"/>	110		Verified	70 - 110
BUN	<input type="checkbox"/>	20		Verified	6 - 20
Creatinine	<input type="checkbox"/>	1.0		Verified	0.5 - 1.2
Sodium Level	<input checked="" type="checkbox"/>	140			136 - 145
Potassium Level	<input type="checkbox"/>	5.0		Verified	3.6 - 5.0
Chloride	<input type="checkbox"/>	101		Verified	101 - 111
CO2	<input type="checkbox"/>	29		Verified	21 - 31
AGAP	<input checked="" type="checkbox"/>	10			4 - 12
Calcium	<input type="checkbox"/>	8.5		Verified	8.5 - 10.5


info: Notice that the AGAP has been recalculated and the flags have been updated.

 **Hit**  **CTRL+A** to add a call comment⁴⁵ to the modified results.⁴⁶

⁴⁵ This is simply to document the call.

⁴⁶ For more information:  [Refer to the COMMENTS Documentation.](#)

WARNING: The Correction comment must be entered **before** proceeding with the next step. Otherwise, the result will need to be Corrected again to add the comment.


 **Click**  when finished.


Procedure		Result	Flags	Status	Reference Ran...
Glucose Level	<input type="checkbox"/>	110		Verified	70 - 110
BUN	<input type="checkbox"/>	20		Verified	6 - 20
Creatinine	<input type="checkbox"/>	1.0		Verified	0.5 - 1.2
Sodium Level	<input type="checkbox"/>	140		Corrected	136 - 145
Potassium Level	<input type="checkbox"/>	5.0		Verified	3.6 - 5.0
Chloride	<input type="checkbox"/>	101		Verified	101 - 111
CO2	<input type="checkbox"/>	29		Verified	21 - 31
AGAP	<input type="checkbox"/>	10		Corrected	4 - 12
Calcium	<input type="checkbox"/>	8.5		Verified	8.5 - 10.5


info: Notice that the STATUS has been changed to CORRECTED.

FINALLY, LEAVE CORRECTION MODE. This step is very important; it will cause frustration if you don't.

 *Click* **MODE** from the menu bar.

 *Click* **CORRECTION**.




THE CORRECTION MODE ICON  should no longer be visible in the lower left corner of ARE.

⚠ IMPORTANT: Don't forget this step. It will prevent you from entering new results.


VIEWING ORIGINAL RESULTS

When results are modified in Cerner, all of the original information is saved.

 *Open* the accession number in Order Result Viewer (ORV).⁴⁷

⁴⁷  Refer to the **ORDER RESULT VIEWER** Documentation.


 *Double-Click* on the order.



Order	Date	Personnel ID	Status	Procedure	Result	C	Normal Values	Units
BMPNL	11/2/2015 10:34	Perry, Thomas J MLS	Verified	Glucose Level	110		70 - 110	mg/dL
BMPNL	11/2/2015 10:34	Perry, Thomas J MLS	Verified	BUN	20		6 - 20	mg/dL
BMPNL	11/2/2015 10:34	Perry, Thomas J MLS	Verified	Creatinine	1.0		0.5 - 1.2	mg/dL
BMPNL	11/2/2015 11:07	Perry, Thomas J MLS	Corrected	Sodium Level	140		136 - 145	mmol/L
BMPNL	11/2/2015 10:34	Perry, Thomas J MLS	Verified	Potassium Level	5.0		3.6 - 5.0	mmol/L
BMPNL	11/2/2015 10:34	Perry, Thomas J MLS	Verified	Chloride	101		101 - 111	mmol/L
BMPNL	11/2/2015 10:34	Perry, Thomas J MLS	Verified	CO2	29		21 - 31	mmol/L
BMPNL	11/2/2015 11:07	Perry, Thomas J MLS	Corrected	AGAP	10		4 - 12	
BMPNL	11/2/2015 10:34	Perry, Thomas J MLS	Verified	Calcium	8.5		8.5 - 10.5	mg/dL

 *Select* the CORRECTED result.

 *Click* 




The 'Result History' window displays the following information:

Accession: 15-306-000014
Order: BMPNL
Procedure: Sodium Level

Status	Date	Personnel ID	Result	Normal Values	Units
Verified	11/2/2015 10:34	Perry , Thomas J II 150 H		136 - 145	mmol/L
Corrected	11/2/2015 11:07	Perry , Thomas J II 140		136 - 145	mmol/L

The window includes a 'Close' button at the bottom right.

THIS WINDOW WILL show all the modifications which have been made to the result.

 **IMPORTANT:** The current result will appear at the bottom of this list. It is sorted from oldest to newest.


INSTRUMENT QUEUE MODE

When Cerner gets results from the instruments that cannot be AUTO-VERIFIED it will send them to the INSTRUMENT QUEUE.

The INSTRUMENT QUEUE mode is used to review questionable results before VERIFYING them.

This mode will be useful in departments where many results are sent from the analyzers.⁴⁸

⁴⁸ Chemistry, Coagulation, and possibly Urinalysis.




OPENING INSTRUMENT QUEUE

Open ARE.


Click **MODE** from the menu bar.

Select **INSTRUMENT QUEUE MODE**.



COSMETIC DIFFERENCES

The difference between INSTRUMENT QUEUE and ACCESSION MODE is the selection box.



TEST SITE: The Instrument Queue you'd like to view.





PROCEDURE: This field is used to limit the results to a specific test.⁴⁹

⁴⁹ This is an optional field.

ACCESSION: This field shows the accession number you're viewing.⁵⁰


⁵⁰ If it's available.

USING INSTRUMENT QUEUE

Enter the  INSTRUMENT,  BENCH,  SECTION or  SUBSECTION.⁵¹

⁵¹ For more information: [Refer to: ROUTING pg. 45](#)

Click 



THE QUEUE WILL remain here until results are sent over which needs to be viewed.

Procedure	Result	Flags	Status	Reference Range	Service Resour...	Prev	Time
Sodium Level	145	D	Performed	136 - 145	SMC DXC 1	135 L	65 minutes
Potassium Level	4.8		Autoverified	3.6 - 5.0	SMC DXC 1	4.8	65 minutes
Chloride	105		Autoverified	101 - 111	SMC DXC 1	105	65 minutes
CO2	24		Autoverified	21 - 31	SMC DXC 1	21	65 minutes
AGAP	16	H	Performed	4 - 12	SMC DXC 1	9	65 minutes

Note: This Screenshot was taken before Autoverify rules were updated.

Review the results.

Enter any appropriate comments.

Click

MOVING TO THE NEXT RESULT

If you're not ready to VERIFY the result you can move to the next item in the queue.

Click the button on the tool bar.

GO TO THE START OF THE QUEUE

When you reach the end of the list, the queue will be blank.

Click then,

Click

THIS WILL BRING you back to the start of the queue.


CANCELING ORDERS

If the results are questionable, we can cancel the orders using ORV or Pending Inquiry.

ADDITIONAL INFORMATION

GETTING ASSAY INFORMATION

Additional information for each assay can be found by **RIGHT-CLICKING** on the RESULT CELL and selecting the appropriate option from the menu.




CALCULATION EQUATIONS

To view the equation used in calculated values:⁵²

⁵² This will not work for calculations that are performed by rules e.g. *Teg Platelet mapping*.

 **Select SHOW EQUATION...**


A POP-UP WINDOW will appear with the equation and a list of the equation's components.



PROCEDURE INFORMATION

To view the various ranges, delta calculations, units of measure and more:

 **Select PROCEDURE INFORMATION...**



The dialog box is titled "Procedure Information" and contains the following fields:

- Procedure:** Short Description: **AGAP**
- Test Site:** SMC Remisol
- Default Result Type:** Calculation
- Current Result Type:** Calculation
- Delta:** Type: (Not Defined), Value: (Not Defined), Minutes: (Not Defined)
- Numeric:** Units of Measure: (Not Applied), Normal Low: 4, Normal High: 12, Critical Low: (Not Applied), Critical High: 20, Review Low: 0, Review High: 20, Linear Low: (Not Applied), Linear High: (Not Applied), Feasible Low: (Not Applied), Feasible High: (Not Applied)
- Alpha:** Normal: (Not Defined), Responses: (None)

A "Close" button is located at the bottom right.

⚠ IMPORTANT: Since many of the ranges can change based on **SEX** and **AGE**, the information displayed is specific to the current patient. The ranges may not apply to other patients.

PREVIOUS RESULTS

The most recent previous results will appear to right of the result.

TIME SINCE LAST RESULT


PREVIOUS RESULT

Procedure	Result	Flags	Status	Reference Range	Service ...	Prev	Time
Glucose Level			Pending	70 - 110			
BUN			Pending	6 - 20			
Creatinine			Pending	0.5 - 1.2		15.0 H	61 days
Sodium Level			Pending	136 - 145		136	47 days
Potassium Level			Pending	3.6 - 5.0		5.0	47 days
Chloride			Pending	101 - 111		101	47 days
CO2			Pending	21 - 31		28	47 days
AGAP			Pending	4 - 12		7	47 days
Calcium			Pending	8.5 - 10.5			


Notice that the TIME column has two different values. This means that the previous results are mostly likely separate orders.

MULTIPLE PREVIOUS


To VIEW ALL the previous results for an assay:

 *Right Click* on the result cell.

 *Click* **PREVIOUS RESULTS...**



THIS WILL OPEN the PREVIOUS RESULTS window for that assay.


A screenshot of the 'Previous Results: Creatinine' window. It displays a table of historical test results for Creatinine. The table has columns: Collect Date/Time, Result, Accession, Perform Date/Time, Verify Date/Time, Order Date/Time, Procedure, and Result Status. The data shows multiple tests over time, with some results marked as 'Corrected' and others as 'Verified'. A 'Close' button is at the bottom right.

Collect Date/Time	Result	Accession	Perform Date/Time	Verify Date/Time	Order Date/Time	Procedure	Result Status
8/19/2015 13:24	15.0 H	15-231-000016	8/19/2015 13:25	8/19/2015 13:25	8/19/2015 13:24	Creatinine	Verified
8/18/2015 13:16	15.0 H	15-230-000036	8/19/2015 13:21	8/19/2015 13:24	8/19/2015 13:21	Creatinine	Verified
8/11/2015 11:56	1.0	15-223-000018	8/11/2015 11:56	8/11/2015 11:56	8/11/2015 11:56	Creatinine	Verified
7/31/2015 12:40	1.0	15-212-000027	7/31/2015 12:44	7/31/2015 12:44	7/31/2015 12:40	Creatinine	Verified
7/31/2015 11:40	25.0 H	15-212-000028	7/31/2015 12:45	7/31/2015 12:45	7/31/2015 12:41	Creatinine	Verified
7/31/2015 11:32	>25.0 HN	15-212-000024	7/31/2015 13:38	7/31/2015 13:38	7/31/2015 11:29	Creatinine	Verified
7/31/2015 9:57	>25.0 HN	15-212-000013	7/31/2015 10:02	7/31/2015 10:02	7/31/2015 9:39	Creatinine	Verified
7/28/2015 9:04	1.2	15-209-000013	7/28/2015 9:06	7/28/2015 9:06	7/28/2015 9:04	Creatinine	Verified
7/27/2015 9:06	1.2	15-208-000058	7/28/2015 9:11	7/28/2015 9:11	7/28/2015 9:11	Creatinine	Verified
7/17/2015 13:54	21.0 H	15-198-000015	7/17/2015 14:01	7/17/2015 14:01	7/17/2015 13:55	Creatinine	Verified
7/9/2015 12:45	15.0 H	15-190-000053	7/9/2015 12:59	7/9/2015 13:46	7/9/2015 12:59	Creatinine	Corrected
6/25/2015 11:12	0.8	15-176-000031	6/25/2015 11:12	6/25/2015 11:13	6/25/2015 11:12	Creatinine	Corrected
6/17/2015 10:02	1.5 H	15-168-000014	6/17/2015 10:57	6/17/2015 10:57	6/17/2015 10:02	Creatinine	Verified
5/18/2015 11:53	5.0 H	15-138-000038	5/18/2015 11:54	5/18/2015 11:54	5/18/2015 11:52	Creatinine	Verified

TROUBLESHOOTING

Too Many Digits

When trying to modify a result, you may get an alert stating that “Too many Digits were entered.”



THIS OFTEN HAPPENS WHEN changing a result value. When you DOUBLE CLICK on a result cell the cursor is placed to the far right.⁵³

⁵³ In the example, it's to the right of ____5.0|

AGAP	<input type="checkbox"/>			Pending
Calcium	<input checked="" type="checkbox"/>	____ 5.0	LCR	Pending

Cerner treats the leading underscores as zeros.


THE EASIEST FIX is to hit  **BACKSPACE** to clear the results.

AGAP	<input type="checkbox"/>			Pending
Calcium	<input checked="" type="checkbox"/>	_____	LCR	Pending

DECIMALS NOT ALLOWED

The number of places after the decimal point is defined for each assay.

If no decimals have been defined, Cerner will prevent them from being entered.



UNABLE TO ENTER RESULTS

There are two reasons why ARE won't let you enter results.

- The sample isn't IN-LAB
- ARE is in Correction Mode.

THE SAMPLE ISN'T IN-LAB

Results can only be entered if status is PENDING. This means that the sample has been LOGGED-IN to the proper location.⁵⁴

IN THIS EXAMPLE the SER HCG QL cannot be resulted.

Procedure		Result	Flags	Status	Reference Range
Sodium Level	<input checked="" type="checkbox"/>	140		Pending	136 - 145
Potassium Level	<input checked="" type="checkbox"/>	4.0		Pending	3.6 - 5.0
Chloride	<input checked="" type="checkbox"/>	108		Pending	101 - 111
CO2	<input checked="" type="checkbox"/>	22		Pending	21 - 31
AGAP	<input checked="" type="checkbox"/>	10		Pending	4 - 12
Ser hCG QL	<input type="checkbox"/>			Not In-Lab	

LYTES PANEL: Has a status of PENDING. This **CAN** be resulted.


SER HCG QL: Has a status of NOT IN-LAB. This **CANNOT** be resulted.

To FIX THIS the sample must be LOGGED-IN. If you currently have the container, this can be done from within ARE.

Click  from the TOOL BAR.

Log-in the sample.⁵⁵

YOU'RE STILL IN CORRECTION MODE


If you can see the  in the bottom left corner of ARE, then you're currently in CORRECTION MODE. New results cannot be entered if ARE is in CORRECTION MODE.


To FIX THIS disable CORRECTION MODE.


Click **MODE** from the menu bar.

⁵⁴ For additional help with this issue: [Refer to the CONTAINER INQUIRY: USING CONTAINER INQUIRY FOR TROUBLESHOOTING Documentation.](#)

⁵⁵ For more information: [Refer to the SPECIMEN LOG-IN PROCEDURE.](#)

 Click **CORRECTION**.



THE CORRECTION MODE ICON  should no longer be visible in the lower left corner of ARE.

ROUTING

WHAT IS ROUTING?

When Accession Numbers are assigned, Cerner will determine a SERVICE RESOURCE⁵⁶ for the containers. The SERVICE RESOURCE is based on the PROCEDURE,⁵⁷ and COLLECTION LOCATION of the container.

⁵⁶ Testing location.

⁵⁷ e. g. CBC, Basic Panel, HbA_{1c}, etc...

THIS GIVES CERNER the ability to update the status of a sample automatically based on the last place it was LOGGED-IN.

FOR EXAMPLE: Let's say a CBC and a HbA_{1c} are sent to SMCA LABORATORY.

When they're LOGGED-IN to SMC LOGIN, the CBC have an IN-LAB status,⁵⁸ while the HbA_{1c} will have PENDING status.⁵⁹



⁵⁸ Meaning, it's available to be tested.

⁵⁹ Meaning, it's been collected but hasn't been sent yet.

IN ADDITION TO TRACKING, the SERVICE RESOURCES is used in Pending Inquiry as the TEST SITE.


HOW DOES IT WORK?



INSTRUMENTS AND BENCHES AND SECTIONS OH MY!

Every test is **ROUTED** to either an  INSTRUMENT or  BENCH.

 INSTRUMENTS are interfaced instruments.⁶⁰

⁶⁰ They can send the results to Cerner.


 BENCHES are where manual tests are performed.

Since there are so many instruments⁶¹ and benches we need a way to organize them. To do this we have  SECTIONS and  SUBSECTIONS.

⁶¹ We should start a band!

 SECTIONS are departments in the lab.⁶²




⁶² Heme, Chem, UA, etc...

 SUBSECTIONS are areas of the sections.⁶³

⁶³ Manual Chemistry, TEGs, Miscellaneous Micro, etc...


ALL THAT TO SAY:

Tests are routed to a  BENCH or  INSTRUMENT.

 BENCHES and  INSTRUMENTS belong to a  SUBSECTION.

 SUBSECTIONS belongs to a  SECTION.



BRACKENRIDGE




LEGEND

- LABORATORY
- SECTION (DEPARTMENT)
- SUB SECTION
- BENCH
- INSTRUMENT


DELL CHILDREN'S




EDGAR B. DAVIS




HAYS




HIGHLAND LAKES






NORTHWEST



SOUTHWEST



WILLIAMSON



SETON MEDICAL CENTER

