

CERNER TRAINING MANUALS

# SPECIMEN LOG-IN



# CONTENTS

INTRODUCTION	1	TIPS AND TRICKS	15
GETTING STARTED	3	Modifying a Single Container . . . . .	15
Setting Location . . . . .	3	Refreshing Data . . . . .	15
Setting Default Printers . . . . .	4	Printing Labels . . . . .	15
LOG-IN BY ACCESSION NUMBER	5	Keyboard Shortcuts . . . . .	16
Understanding the Spreadsheet . . . . .	5	TROUBLESHOOTING	17
Log-in Accessions Procedure . . . . .	6	Previously Entered . . . . .	17
LOG-IN BY LIST	11	Invalid Collection Time . . . . .	17
List Pop-Up Alerts . . . . .	12	Separate Collections . . . . .	18
Recommended Procedure . . . . .	13	I've Logged a Sample into the Wrong Location	18





# INTRODUCTION

Specimen Log-in is the application used to update sample locations within Cerner. This application will most commonly be used to log-in samples which have been collected by the floors.

Specimen Log-in is also used to log-in “Transfer Lists” sent by other facilities.

Cerner will not allow orders to be resulted until they have been LOGGED-IN to the testing laboratory. In most cases, the instruments will not know which tests to perform.


OPEN Specimen Log-in by clicking the  icon from the App-Bar.



THERE ARE THREE way's samples can be logged-in:

**LIST:** Log-in an entire list. Either a *Collection List*, or a *Transfer List*

**PATIENT:** This is to log samples in using a Patient Identifier. In *most* cases, this method is not recommended,<sup>1</sup> so it will not be discussed in this procedure.

**ACCESSION:** Log-in samples using the Accession Number.<sup>2</sup>  Refer to: **LOG-IN BY ACCESSION NUMBER** pg. 5.

<sup>1</sup> However, it may be useful for sites where the Laboratory collects samples.


<sup>2</sup> i. e. *Scanning the bar-code*.



## GETTING STARTED

The first time opening Specimen Log-in, it's a good idea to modify the default values.

 *Open* Specimen Log-in by clicking the  icon from the App-Bar.

 *Select* "Accession" from the "Log-in By" options.


 *Click*  to open "Log-In By ACCESSION."

## SETTING LOCATION

The most important setting in Specimen Log-in is the LOCATION. In most cases, the location will be set to your laboratory's LOG-IN LOCATION.

**Fun Fact:** Most LOG-IN related mistakes come from this being set to the wrong location.

LOCATION



- BH LOGIN
- DC LOGIN
- SEBD LOGIN
- SHC LOGIN
- SHL LOGIN
- SMC LOGIN
- SNW LOGIN
- SSW LOGIN
- SWC LOGIN

*e. g. Dell's log-in location is "DC Login."*


Of course, there are exceptions. Some departments will require an alternate Log-in Location e.g. *Microbiology, Blood Bank, etc....* Departments which require special log-in locations should provide that information.

**⚠ IMPORTANT:** This setting is “Sticky” and there is no “Default” value. If you temporarily change it from your normal log-in location,<sup>3</sup> it needs to be changed back.<sup>4</sup>

<sup>3</sup> e.g. *Receiving microbiology samples.*

<sup>4</sup> This will be a common mistake for people starting out.


 **Select** your LOG-IN LOCATION from the drop-down menu.





## SETTING DEFAULT PRINTERS


 **Click** TASK from the menu-bar.

 **Select** LABEL PRINTER...



 **Set** both drop-down menus to the nearest label printer.

 **Click**  to save the settings.



**Note:** To save the settings temporarily, ☐ **Un-check** the SET AS DEFAULT... options.



# LOG-IN BY ACCESSION NUMBER

Log-in by Accession Number is used to log-in samples using the *barcoded* label on the specimen container.

## UNDERSTANDING THE SPREADSHEET

The SPECIMEN LOG-IN SPREADSHEET lists the containers which have been scanned to log-in.

C...	Accession	Container	Order	Coll Date	Coll Time	Coll ID	Order	Priority	Status	Coll Method	Rec D
<input checked="" type="checkbox"/>	15-268-000019	A : 3 mL Lavender	CBC	10/1/2015			CBC	RT	Dispatch	Venous Draw	10/1/2
<input checked="" type="checkbox"/>		D : 5 mL Blue									
<input checked="" type="checkbox"/>			D-Dimer	10/1/2015			D-Dimer	RT	Dispatch	Venous Draw	10/1/2
<input checked="" type="checkbox"/>			PTT	10/1/2015			PTT	RT	Dispatch	Venous Draw	10/1/2
<input checked="" type="checkbox"/>			PROTIM	10/1/2015			PROTIM	RT	Dispatch	Venous Draw	10/1/2
<input checked="" type="checkbox"/>			FIB	10/1/2015			FIB	RT	Dispatch	Venous Draw	10/1/2
<input checked="" type="checkbox"/>	15-274-000041	A : 100 mL Urine Cup	UAS	10/1/2015	18:29		UAS	RT	Collected	Bag Collection	10/1/2

In this example, three containers have been scanned. Container **A** and **D** of one accession number and container **A** of another.

THE ROWS are grouped by accession numbers and Containers.

As explained in [ACCESSION NUMBERS](#), accession numbers can have multiple containers.

UNDER EACH ACCESSION NUMBER are the containers which have been scanned.<sup>5</sup>

UNDER EACH CONTAINER are the orderables associated with it.

C...	Accession	Container	Order
<input checked="" type="checkbox"/>	15-268-000019		
<input checked="" type="checkbox"/>		A : 3 mL Lavender	
<input checked="" type="checkbox"/>			CBC
<input checked="" type="checkbox"/>		D : 5 mL Blue	
<input checked="" type="checkbox"/>			D-Dimer
<input checked="" type="checkbox"/>			PTT
<input checked="" type="checkbox"/>			PROTIM
<input checked="" type="checkbox"/>			FIB
<input checked="" type="checkbox"/>	15-274-000041		
<input checked="" type="checkbox"/>		A : 100 mL Urine Cup	
<input checked="" type="checkbox"/>			UAS

<sup>5</sup> In the example above we can see that for accession number 15-268-000019, container **A**, and **D** have been scanned.

This image shows the two containers, and their orders.

## CHECK-BOXES

Each row contains a check-box. By default, they are ☒d. These can be used to pick which containers are LOGGED-IN.

**SELECTED** ☒: The orderable or container is selected for an action.<sup>6</sup>

<sup>6</sup> e. g. Log-in, Miss...The two main functions.


**DE-SELECTED** ☐: The orderable or container is not select for an action.


**DISABLED** ☐: The orderable or container has been LOGGED-IN.

AFTER A CONTAINER has been LOGGED-IN, the check-box will become DISABLED and it can no longer be modified.


[Refer to: REFRESHING DATA pg. 15](#) for information on clearing these containers.

## LOG-IN ACCESSIONS PROCEDURE

OPEN Specimen Log-in by clicking the  icon from the App-Bar.

 *Select* ACCESSION from the LOG-IN BY options.


 *Click*  to open LOG-IN BY ACCESSION.


 *Check* that your LOG-IN location is set appropriately.

 Click on the ACCESSION entry field.<sup>7</sup>



<sup>7</sup> Make sure you see that blinky line thing: |  
**info:** Skipping this step after changing locations will cause a minor inconvenience.


 Scan the bar-code of the container.



**! IMPORTANT:** While, the accession number can be manually entered, it's not recommended. Accidentally omitting the container identifier will cause uncollected containers to be received in error.

Besides, scanning the tubes is much easier and faster.


SPECIMEN LOG-IN WILL update with the container scanned.



C...	Accession	Container	Order	Coll Date	Coll Time	Coll ID	Order	Priority	Status	Coll Method
<input checked="" type="checkbox"/>	15-268-000019	A: 3 mL Lavender	CBC	10/1/2015			CBC	RT	Dispatch	Venous Draw

**! IMPORTANT:** If the accession number has multiple containers,<sup>8</sup> scan them as well. Otherwise, only the scanned container will be logged-in.

<sup>8</sup> Since a Green and Lav have different Container Identifiers, they are treated individually.



C...	Accession	Container	Order	Coll Date	Coll Time	Coll ID	Order	Priority	Status	Coll Method	Rec Date	R
<input checked="" type="checkbox"/>	15-268-000019	A: 3 mL Lavender	CBC	10/8/2015			CBC	RT	Dispatch	Venous Draw	10/8/2015	
<input checked="" type="checkbox"/>		B: 3 mL Mint Green	BMPNL	10/8/2015			BMPNL	RT	Dispatch	Venous Draw	10/8/2015	

## UPDATING COLLECT INFORMATION


**⚠ IMPORTANT:** This information should be entered at the time of collection. If it hasn't, the collecting personnel should be notified and instructed to enter it.

IF NEEDED, the collection information can be updated. This should only be done for outreach samples or if there are extenuating circumstances.


**info:** The following instructions will update all the containers.<sup>9</sup> If only one value needs to be update [Refer to: MODIFYING A SINGLE CONTAINER pg. 15.](#)

<sup>9</sup> Except the ones which have already been LOGGED-IN. Their check boxes will be disabled ☐.


 **Click** the COLL ID COLUMN HEADER.




C...	Accession	Container	Order	Coll Date	Coll Time	Coll ID	Order
<input checked="" type="checkbox"/>	15-268-000019						
<input checked="" type="checkbox"/>		A : 3 mL Lavender					
<input checked="" type="checkbox"/>			CBC	10/8/2015			CBC
<input checked="" type="checkbox"/>		B : 3 mL Mint Green					
<input checked="" type="checkbox"/>			BMPNL	10/8/2015			BMPNL

 **Enter** the USER-NAME<sup>10</sup> of the person who collected the sample.

<sup>10</sup> This needs to be the Cerner Username. Generic names exist for outside facilities.



 **Click** the COLL TIME COLUMN HEADER.



C...	Accession	Container	Order	Coll Date	Coll Time	Coll ID	Order
<input checked="" type="checkbox"/>	15-268-000019						
<input checked="" type="checkbox"/>		A : 3 mL Lavender					
<input checked="" type="checkbox"/>			CBC	10/9/2015		TPERR005	CBC
<input checked="" type="checkbox"/>		B : 3 mL Mint Green					
<input checked="" type="checkbox"/>			BMPNL	10/9/2015		TPERR005	BMPNL

 **Enter** the COLLECTION TIME using a 24h format.<sup>11</sup>

<sup>11</sup> This needs to be earlier than the RECEIVE TIME which defaults to Now.

Coll Time

21:21

OKCancel

Click the COLL DATE COLUMN HEADER.

C...	Accession	Container	Order	Coll Date	Coll Time	Coll ID	Order
<input checked="" type="checkbox"/>	15-268-000019						
<input checked="" type="checkbox"/>		A : 3 mL Lavender					
<input checked="" type="checkbox"/>			CBC	10/10/2015	21:21	TPERR005	CBC
<input checked="" type="checkbox"/>		B : 3 mL Mint Green					
<input checked="" type="checkbox"/>			BMPNL	10/10/2015	21:21	TPERR005	BMPNL


Enter the COLLECTION DATE using a mm/dd/yyyy format.<sup>12</sup>

Coll Date


10/09/2015



OKCancel

If Needed, Hit  Ctrl+B to print the labels.

Click  to Log-in<sup>13</sup>

C...	Accession	Container	Order	Coll Date	Coll Time	Coll ID	Order	Priority	Status
<input type="checkbox"/>	15-268-000019								
<input type="checkbox"/>		A : 3 mL Lavender							
<input type="checkbox"/>			CBC	10/9/2015	21:21	TPERR005	CBC	RT	Collected
<input type="checkbox"/>		B : 3 mL Mint Green							
<input type="checkbox"/>			BMPNL	10/9/2015	21:21	TPERR005	BMPNL	RT	Collected


<sup>12</sup>  **IMPORTANT:** This date defaults to TODAY. If the sample was collected around midnight, make sure this date is correct.

<sup>13</sup>  **Hotkey:** CTRL+L will also work.  
**info:** If you get a pop-up asking about “Separate Collections”  Refer to: **SEPARATE COLLECTIONS** pg. 18.

AT THIS POINT, the samples have been LOGGED-IN. Those containers will remain in the window until it is refreshed.

THE MISS BUTTON

Please ignore this button. As of right now it doesn’t fit into our work flow.

 Refer to: **REFRESHING DATA** pg. 15 for information on clearing these containers.



## LOG-IN BY LIST


Samples sent from other facilities need to be LOGGED-IN to the testing laboratory before they can be resultated.


Log-in by List is used to LOG-IN all of the samples contained on a COLLECTION or TRANSFER LIST. This section discusses that process. However, there is a better way to LOG-IN entire lists at the end of this section.


**info:** The recommended way to Log-In lists is by using LOG-IN BY ACCESSION, then use LOG-IN BY LIST to confirm that all the samples have arrived.


 Refer to: [RECOMMENDED PROCEDURE pg. 13](#) for more information.

This will save a considerable amount of time, while also providing a quick and accurate way of verifying that all the samples have arrived.

OPEN Specimen Log-in by clicking the  icon from the App-Bar.




 Click the LIST option.




 Enter the COLLECTION or TRANSFER LIST NUMBER.<sup>14</sup>

<sup>14</sup> Using the appropriate text box.

 **Update** the LIST DATE.<sup>15</sup>

<sup>15</sup> This is the date that the list was created.



PathNet Collections: Specimen Log-In

Task Edit View Help

Log-In By:

☐ User

☐ Location

☒ List

☐ Patient

☐ Accession

Collection List Number:

Transfer List Number:

51

List Date:


10/16/2015


Retrieve Close

Ready A3291 TPERR005 10/16/2015 15:46

 **Click** 

THE LOG-IN BY TRANSFER LIST window will open and the containers of the list will populate the Spreadsheet.<sup>16</sup>

<sup>16</sup>  Refer to: UNDERSTANDING THE SPREAD-SHEET pg. 5 for more information.



Lab 32 @ 10/16/2015 - From SMC Login To To BH - Log-In By TRANSFER LIST

Task View Help

Demographics

Name: ZZTEST, UMCBSCEM6

MRN: 5026159

Loc: Brackenridge Hospital / Discharged

FIN: 5503248150

Gender: Female

Ordered by:

DOB: 1/1/1984


Admitted: 4/8/2014 10:46


C...	Accession	Container	Order	Coll Date	Coll Time	Coll ID	Order	Priority	Status	Coll Method	Rec Date	Rec Time	Cont/Vol
	15-270-000115	A : 3 mL Mint Green	CMPNL	9/27/2015	8:52		CMPNL	RT	Collected	Venous Draw	10/16/2015	18:07	3 mL Mint G
	15-147-000022	A1FP	A1FP	9/27/2015	8:52		A1FP	RT	Collected	Venous Draw	10/16/2015	18:07	3 mL Mint G
	15-078-000026	A : 100 mL Urine Cup	YUPEP	5/27/2015	12:27		YUPEP	RT	Collected	Container/Tube	10/16/2015	18:07	100 mL Urine
	15-147-000024	G : 5 mL Gold	RUBE	3/19/2015	15:41		RUBE	RT	Collected	Venous Draw	10/16/2015	18:07	5 mL Gold
	15-268-000035	B : 5 mL Gold	RUBE	5/27/2015	12:33		RUBE	RT	Collected	Venous Draw	10/16/2015	18:07	5 mL Gold
	15-268-000035	B : 5 mL Gold	A1FP	9/25/2015	16:36		A1FP	RT	Collected	Venous Draw	10/16/2015	18:07	5 mL Gold



Location: BH Login


Log In Miss Details



Ready A3291 TPERR005 10/16/2015 18:07

 **Compare** the list of containers on the Spreadsheet with the containers received.<sup>17</sup>

<sup>17</sup> This process is time consuming.  Refer to: RECOMMENDED PROCEDURE pg. 13 for a quicker, and more accurate method.


 **Check**  the containers received.<sup>18</sup>

<sup>18</sup> Don't  containers if they have not been sent. These samples are considered "Missing."

 **Click**  to Log-in.<sup>19</sup>

<sup>19</sup>  Hotkey: CTRL+L will also work.


## LIST POP-UP ALERTS

**info:** If you get a pop-up asking about "Separate Collections"  Refer to: SEPARATE COLLECTIONS pg. 18.

LOG-IN BY LIST has two Pop-Up boxes it will display.



LIST NOT FOUND will be displayed if it cannot find a list with that number, for that date.



WHAT IF THE LIST EXISTS? There are three things that can cause this Pop-Up:


**WRONG TYPE:** The LIST NUMBER was entered into the wrong field.<sup>20</sup>

<sup>20</sup> e.g. A Transfer List Number was entered into the Collections List field.

**WRONG DATE:** The DATE FIELD does not have the correct date.

**WRONG NUMBER:** The LIST NUMBER was not entered correctly.

NO ORDERS FOUND will be displayed if all of the containers for a list have been LOGGED-IN.





THIS MEANS all the samples have been LOGGED-IN. The list is now "Empty."

## RECOMMENDED PROCEDURE


This method eliminates the need to manually compare the TRANSFER LIST with the containers received.

It will save a considerable amount of time, and as long as the barcode scanners are used, it will identify any missing samples.

 *Log-in* each container using **LOG-IN BY ACCESSION**.


 *Attempt* to LOG-IN the LIST using **LOG-IN BY LIST**.

IF THERE ARE NO MISSING samples, LOG-IN BY LIST will display an alert.



**⚠ WARNING:** *This alert looks dangerously close to the **LIST NOT FOUND** alert.*

IF SAMPLES ARE MISSING, the LOG-IN BY TRANSFER LIST window will open and the missing containers will populate the Spreadsheet.



Log-In: 52 @ 10/16/2015 - From SMC Login To To BH - Log-In By TRANSFER LIST

Task View Help

Demographics  
 Name: 221TEST, UMCSCEN6  
 FIN: 5503248150  
 DOB: 1/1/1984  
 MRN: 5026159  
 Gender: Female  
 Loc: Brackenridge Hospital / Discharged  
 Ordered by:  
 Admitted: 4/8/2014 10:45

C	Accession	Container	Order	Col Date	Col Time	Col ID	Order	Priority	Status	Col Method	Rec Date	Rec Time	Cont/Vol
	15-270-000015	A : 3 mL Mint Green											
			CMPNL	9/27/2015	8:52		CMPNL	RT	Collected	Venous Draw	10/16/2015	18:07	3 mL Mint G
			A1FP	9/27/2015	8:52		A1FP	RT	Collected	Venous Draw	10/16/2015	18:07	3 mL Mint G
	15-147-000022	A : 100 mL Urine Cup											
			YUPEP	5/27/2015	12:27		YUPEP	RT	Collected	Container/Tube	10/16/2015	18:07	100 mL Urine
	15-078-000026	G : 5 mL Gold											
			RUBE	3/19/2015	15:41		RUBE	RT	Collected	Venous Draw	10/16/2015	18:07	5 mL Gold
	15-147-000024	B : 5 mL Gold											
			RUBE	5/27/2015	12:33		RUBE	RT	Collected	Venous Draw	10/16/2015	18:07	5 mL Gold
	15-268-000035	B : 5 mL Gold											
			A1FP	9/25/2015	16:36		A1FP	RT	Collected	Venous Draw	10/16/2015	18:07	5 mL Gold

Location: BH Login Log In Miss Details

Ready A3291 TPERR005 10/16/2015 18:07

These samples may have been missed during the initial LOG-IN process, or they may still be at their site of origin. Either way, it's time for some detective work.

# TIPS AND TRICKS

Here are some Specimen Log-in tips. Knowing these functions can save time, and simplify your work-flow.

## MODIFYING A SINGLE CONTAINER

If only one value needs to be modified, it can be done by clicking on its cell in the SPECIMEN LOG-IN SPREADSHEET.

C...	Accession	Container	Order	Coll Date	Coll Time	Coll ID	Order
<input checked="" type="checkbox"/>	15-268-000019	D : 5 mL Blue					
<input checked="" type="checkbox"/>			D-Dimer	10/15/2015			D-Dimer
<input checked="" type="checkbox"/>			PTT	10/15/2015			PTT
<input checked="" type="checkbox"/>			PROTIM	10/15/2015			PROTIM
<input checked="" type="checkbox"/>			FIB	10/15/2015			FIB

The downside to this method is the collection information needs to be updated for each orderable. This can be time consuming.

## REFRESHING DATA

After the samples have been LOGGED-IN, they will remain in the SPECIMEN LOG-IN SPREADSHEET. While they don't need to be removed before continuing with other samples, it is a good idea to periodically clear them out.

Click the REFRESH DATA<sup>21</sup>  icon from the toolbar.

<sup>21</sup> The glasses and cup.



**▲ IMPORTANT:** This will only remove samples which have been **Logged-in**.<sup>22</sup> Unfortunately, you can't remove samples that have not been **Logged-in** or **Missed**<sup>23</sup> without closing the window and re-opening it.

<sup>22</sup> Their Check-boxes will be disabled.

<sup>23</sup> Sorry.


## PRINTING LABELS

Before samples have been LOGGED-IN, the labels can be quickly reprinted.

Hit  **CTRL+B**

KEYBOARD SHORTCUTS

SPECIMEN LOGIN SHORTCUTS	
 CTRL+L	LOG IN SELECTED SPECIMEN
 CTRL+B	PRINT LABELS


Table 1:  Specimen Login Keyboard Shortcuts

# TROUBLESHOOTING

## PREVIOUSLY ENTERED

If a container is currently in the SPECIMEN LOG-IN SPREADSHEET it cannot be scanned in.

This is one of the most frustrating things about Specimen Log-in.




To FIX this, the SPREADSHEET needs to be refreshed. [Refer to: REFRESHING DATA pg. 15](#) for information on how that is done.


## INVALID COLLECTION TIME

Cerner will not accept COLLECTION TIMES that are more recent than the RECEIVED TIME.

Sorry time travelers.



To FIX this, simply:

 *Correct* the times to their actual values to continue.

## SEPARATE COLLECTIONS

Occasionally, containers will be collected at different times.<sup>24</sup>

Since the containers have their own collection information, this is acceptable. Cerner will display a *friendly* alert asking if this is correct.

<sup>24</sup> e. g. Blue top was collected at 15:00, Green top was collected at 16:00.

Entry	Accession	Container	Coll Date	Coll Time	Coll ID	Coll
Your Entry	15-268-000019	C : 3 mL Lavender	10/15/2015	18:59	TPERR005	Veno
Previous Log-In	15-268-000019	A : 3 mL Lavender	10/9/2015	21:21	TPERR005	Veno
	15-268-000019	B : 3 mL Mint Green	10/9/2015	21:21	TPERR005	Veno

Were these really separate collections?

Yes No Skip

IF THE INFORMATION is correct:<sup>25</sup>

Click the **Yes** button.

<sup>25</sup> The samples were collected at two separate times.

IF THE SAMPLES were collected at the same time:

Click the **No** button.

**⚠ WARNING:** This option will apply the collection information from the first container to the container you're trying to Log-in. Only do this if you've entered the wrong collection time in the previous steps.

## I'VE LOGGED A SAMPLE INTO THE WRONG LOCATION

DON'T PANIC! Simply log it into the correct location.

