

CERNER TRAINING MANUALS

# TRANSFER SPECIMENS



# CONTENTS


INTRODUCTION	1	UPDATING TRANSFER LISTS	9
		Adding or Removing	
GETTING STARTED	3	Containers . . . .	9
Status . . . . .	3	Deleting a Transfer List .	10
Route . . . . .	4	Transferring Transfer Lists	11
CREATE A TRANSFER LIST	5	TROUBLESHOOTING	13
List Build . . . . .	5	I Can't Find a Container	13
		I Can't Find a List . . . .	15





# INTRODUCTION

When testing is performed at another facility, it needs to be TRANSFERED in Cerner before it's picked up by the carrier.

This process tells Cerner that the sample is on its way from the one location to the next. It will update the sample status in Container Inquiry and help with specimen tracking.

OPEN TRANSFER SPECIMENS by clicking the  icon from the App-Bar.<sup>1</sup>

<sup>1</sup>  Refer to the APP-BAR PROCEDURE if you need help adding it.



PathNet Collections: Transfer Specimens

Task View Help

Transfer Lists

☒ Not transferred ☐ Transferred

Locations and Date

From: To: Date: Retrieve

From SMC Login To BH 11/04/2015

Lists

Accession	Order	Patient Name	Patient ID	Desc	Seq

New... Modify... Delete Transfer


No Transfer Lists Found A3291 TPERR005 11.28

TRANSFER LISTS can be thought of as a log of samples going from one location to another.



# GETTING STARTED

This section will discuss the main window and how to read it.



**STATUS:** changes which TRANSFER LISTS are viewed.<sup>2</sup>

<sup>2</sup> TRANSFERED VS NOT TRANSFERED

**ROUTE:** is used to define the SOURCE, DESTINATION and DATE.

**LISTS:** shows the lists meeting the criteria set in **STATUS** and **ROUTE**.

**ITEMS:** shows the CONTAINERS on the selected LIST.

## Status

When a TRANSFER LIST is created, its state will be NOT TRANSFERED. The list can be modified or deleted.

When the samples are picked up by the courier, the list needs to be TRANSFERED. This will put it in a TRANSFERED state.

Finally, when the samples have been LOGGED-IN to the destination the TRANSFER LIST will be cleared.

Table 2: Transfer List States

STATES	WHAT IT MEANS
NOT TRANSFERRED	The samples are still at the original lab.
TRANSFERRED	The sample is in route to the destination.
CLEARED	The samples have arrived at the destination.

## Route

When viewing TRANSFER LISTS, these options are used to define which lists are displayed.

**FROM:** Your current location.<sup>3</sup>

**To:** The destination.<sup>4</sup>

**DATE:** The date the TRANSFER LIST was CREATED.

**▲ IMPORTANT:** THE **FROM** LOCATION IS WHERE YOU ARE SENDING THE SAMPLE.  
*e.g. If the test is performed by Arup, but needs to go to BH REFERRALS first, use **BH REFERRALS** as the DESTINATION.*

<sup>3</sup> Always choose the option beginning with **FROM**.

<sup>4</sup> *e.g.* To BH, To BH MICRO, To BH REFERRALS, etc...



# CREATE A TRANSFER LIST


There are two ways to build TRANSFER LISTS: LIST BUILD and MANUAL BUILD.


LIST BUILD will allow you to choose from the containers which need to be sent.

MANUAL BUILD allows you to build a list by scanning the containers.<sup>5</sup>

<sup>5</sup> This method is the best way to make sure the list is accurate.


## List Build

OPEN TRANSFER SPECIMENS by clicking the  icon from the App-Bar.<sup>6</sup>

<sup>6</sup>  Refer to the APP-BAR PROCEDURE if you need help adding it.

 Click <sup>7</sup>

<sup>7</sup> Bottom right corner of the screen.




 *Select* LIST BUILD mode.

Mode

☒ List build ☐ Manual build

 *Select* the **FROM LOCATION** from the following options:



- FROM BERTRAM                      - FROM LOCKHART                      - FROM SHOAL CRK
- FROM BH LOGIN                      - FROM LOCKHART S                      - FROM SMC LOGIN
- FROM BURNET OP                      - FROM MARBLE FAL                      - FROM SNW LOGIN
- FROM DC LOGIN                      - FROM MCCARTHY                      - FROM SP CHILDR
- FROM KINGSLAND                      - FROM SEBD LOGIN                      - FROM SSW LOGIN
- FROM KOZMETSKY                      - FROM SHC LOGIN                      - FROM SWC LOGIN
- FROM LAMPASAS                      - FROM SHL LOGIN                      - FROM TOPFER

 *Select* the **TO LOCATION**.<sup>8</sup>

**⚠ IMPORTANT:** THE SELECTION MUST BEGIN WITH THE WORD **To**. THE MOST COMMON CHOICES WILL BE: To BH, To BH REFERRALS AND To BH MICRO.<sup>9</sup>

<sup>8</sup> This is the location the carrier will be delivering the samples.

<sup>9</sup> The other options are there for the referrals department.


 *Click* 


From Available				
	Accession	Order	Patient Name	Patient ID
1	15-264-000058A	BUNU	TESTING, NOTWC	1017792
2	15-293-000039A	HBA1C	TEST, VON TESTINGSTON	603834200
3	15-293-000042D	RPR	TEST, VON TESTINGSTON	603834200
4	15-308-000015A	HBA1C	TEST, VON TESTINGSTON	603834200
5	15-308-000018A	TSH	TEST, VON TESTINGSTON	603834200
6	15-309-000033A	CEA	ZZTEST, SUT15AOREG1	1017816

This pulls up a list of samples which are available to be sent to the destination.

**info:** The TSH is listed here because it was manually **RE-ROUTED** to Brackenridge using PENDING INQUIRY.



 *Select* the containers to be added to the list.<sup>10</sup>


<sup>10</sup> Pressing  **CTRL** while **CLICKING** allows you to select multiple containers.

 *Click*  to add the containers to the new list.<sup>11</sup>

<sup>11</sup> The samples will move to the new list, and are no longer available to be put in another list.

To Selected				
	Accession	Order	Patient Name	Patient ID
1	15-293-000042D	RPR	TEST, VON TESTINGSTON	6038342000 3 ml
2	15-308-000015A	HBA1C	TEST, VON TESTINGSTON	6038342000 3 ml
3	15-308-000018A	TSH	TEST, VON TESTINGSTON	6038342000 3 ml

 Click  when finished.



 Click 

THE NEW LIST now appears in the LISTS AREA and its containers can be viewed when clicked.


Lists		From:		To:		Transfer Temp:	
2		From SMC Login		To BH		(All)	
6		Accession	Order	Patient Name	Patient ID	Desc	Seq
		1	15-293-000042D RPR	TEST, VON TESTINGSTON	6038342000	3 mL SST	1
		2	15-308-000015A HBA1C	TEST, VON TESTINGSTON	6038342000	3 mL Lavender	2
		3	15-308-000018A TSH	TEST, VON TESTINGSTON	6038342000	3 mL Mint Green	3



# UPDATING TRANSFER LISTS

Transfer lists which have not been TRANSFERED can be modified. This means containers can be added or removed from the list.

## Adding or Removing Containers


 *Select* **NOT TRANSFERED** from the TRANSFER LIST options.




Transfer Lists


☒ Not transferred ☐ Transferred

 *Select* the **FROM LOCATION**.

 *Select* the **TO LOCATION**.

 *Update* the **DATE**.

 *Click* 



PathNet Collections: Transfer Specimens

Task View Help

Transfer Lists

☒ Not transferred ☐ Transferred


Locations and Date

From:  To:  Date:

Lists

	From: From SMC Login	To: To BH	Transfer Temp: (All)
2	Accession	Order	Patient Name
6	15-293-000042D	RPR	TEST, VON TESTINGSTON
	15-308-000015A	HBA1C	TEST, VON TESTINGSTON
	15-308-000018A	TSH	TEST, VON TESTINGSTON
			Patient ID
			6038342000
			Desc
			3 mL SST
			Seq
			1
			2
			3


Ready A3291 TPERR005 14:11

 *Select* the list to modify.

Click 

Move any containers to or from the list.

**info:** Removing an item from a transfer list will make it available to be put on another list.

Click  when finished.

## Deleting a Transfer List

Select **NOT TRANSFERED** from the TRANSFER LIST options.

Select the **FROM LOCATION**.

Select the **TO LOCATION**.

Update the **DATE**.

Click 

Accession	Order	Patient Name	Patient ID	Desc	Seq
15-293-000420	RPR	TEST, VON TESTINGSTON	6038342000	3 mL SST	1
15-308-000015A	HBAIC	TEST, VON TESTINGSTON	6038342000	3 mL Lavender	2
15-308-000018A	TSH	TEST, VON TESTINGSTON	6038342000	3 mL Mint Green	3

*Select* the list to delete.

*Click*

**⚠ WARNING:** THE TRANSFER LIST WILL INSTANTLY BE DELETED, BUT FRET NOT, IT CAN EASILY BE RECREATE IT IF NEEDED.

## Transferring Transfer Lists

*Select* **NOT TRANSFERED** from the TRANSFER LIST options.

*Select* the **FROM LOCATION**.

*Select* the **TO LOCATION**.

*Update* the **DATE**.

*Click*

PathNet Collections: Transfer Specimens

Task View Help

Transfer Lists

☒ Not transferred ☐ Transferred

Locations and Date

From: From SMC Login To: To BH Date: 11/05/2015 Retrieve

Lists

From:	To:	Transfer Temp:
From SMC Login	To BH	(AB)
Accession	Order	Patient Name
1	15-293-0000420 RPR	TEST_VON TESTINGSTON
2	15-308-000015A HBAC	TEST_VON TESTINGSTON
3	15-308-000018A TSH	TEST_VON TESTINGSTON

Ready A3291 TPERR005 14:11

Select the list to transfer.

Click

### View Transferred Lists

Select **TRANSFERED** from the TRANSFER LIST options.

PathNet Collections: Transfer Specimens

Task View Help

Transfer Lists

☐ Not transferred ☒ Transferred

Select the **FROM LOCATION**.

Select the **TO LOCATION**.

Update the **DATE**.

Click

PathNet Collections: Transfer Specimens

Task View Help

Transfer Lists

☐ Not transferred ☒ Transferred

Locations and Date

From: From SMC Login To: To BH Date: 11/05/2015 Retrieve

Lists

From:	To:	Transfer Temp:
From SMC Login	To BH	(AB)
Accession	Order	Patient Name
1	15-293-0000420 RPR	TEST_VON TESTINGSTON
2	15-308-000015A HBAC	TEST_VON TESTINGSTON
3	15-308-000018A TSH	TEST_VON TESTINGSTON

Ready A3291 TPERR005 14:44



# TROUBLESHOOTING

TRANSFER SPECIMENS can be a confusing application. This section will break down into two groups: **I CAN’T FIND A CONTAINER** and **I CAN’T FIND A LIST**.

## *I Can’t Find a Container*

Three things can prevent a container from appearing in the AVAILABLE CONTAINERS LIST.


From Available					
	Accession	Order	Patient Name	Patient ID	
1	15-264-000058A	BUNU	TESTING, NOTWC	1017792	100
2	15-293-000039A	HBA1C	TEST, VON TESTINGSTON	6038342000	3 ml


## *Check the List Options*

Many times, this issue is caused because either the **To LOCATION, FROM LOCATION** or **DATE** are incorrectly entered.

Sometimes, it can be caused because a **TRANSPORT TEMPERATURE** is set incorrectly.<sup>12</sup>

<sup>12</sup> Setting this to all is a way to make sure that’s not the issue.


 *Check* the **FROM** Location.

 *Check* the **To** Location.

 *Select* **(ALL)** in the **TRANSFER TEMP:** dropdown.


 *Select* **ALL CONTAINERS**


 *Click* 

Locations		To:		Transfer Temp:		<input checked="" type="radio"/> All containers <input type="radio"/> Aliquots only <input type="radio"/> Specimens only			
From: From SMC Login		To: To BH		(All)					
From Available									
	Accession	Order	Patient Name	Patient ID					
1	15-264-000058A	BUNU	TESTING, NOTWC	1017792	100				
2	15-293-000039A	HBA1C	TEST, VON TESTINGSTON	6038342000	3 ml				
						To Selected			
	Accession	Order	Patient Name						

## Check Container Status

 *Open*  Container Inquiry.

 *Enter* the accession number in the ACCESSION field.

 *Check* the status of the sample.<sup>13</sup>

<sup>13</sup> Is the sample received at your laboratory?

STATUS LOCATION

Event	Date/Time	User	Location	Detail	Update Date/Time	Update User
Collected	11/4/2015 11:25	(0)	SMC Login		11/4/2015 11:25	SYSTEM
Received	11/4/2015 11:25	TPERR005	SMC Login		11/4/2015 11:25	SYSTEM
Logged Out	11/5/2015 14:39	TPERR005	SMC Login		11/5/2015 14:39	TPERR005
In Transit	11/5/2015 14:39	TPERR005			11/5/2015 14:39	TPERR005

THE MOST RECENT status should be **RECEIVED**, and the location should be your **LOGIN LOCATION**.

SAMPLE ALREADY SENT:<sup>14</sup>

<sup>14</sup> This tells us that the sample already placed on a TRANSFER LIST.  
Nothing needs to be done.

Event	Date/Time	User	Location	Detail
Collected	11/4/2015 11:25	(0)	SMC Login	
Received	11/4/2015 11:25	TPERR005	SMC Login	
Logged Out	11/5/2015 14:39	TPERR005	SMC Login	
In Transit	11/5/2015 14:39	TPERR005		

SAMPLE HASN'T BEEN RECEIVED:<sup>15</sup>

<sup>15</sup> The last STATUS is not "Received" at "SMC Login."  
This sample needs to be LOGGED-IN using Specimen Log-in.

Event	Date/Time	User	Location
Dispatched	11/4/2015 6:01	SYSTEM	

THIS IS A **GOOD** STATUS:<sup>16</sup>

<sup>16</sup> We may need to check the **ROUTING**.

Event	Date/Time	User	Location	Detail
Collected	11/4/2015 11:25	(0)	SMC Login	
Received	11/4/2015 11:25	TPERR005	SMC Login	


### Check the Routing

Open  Container Inquiry.

Enter the accession number in the ACCESSION field.

Right Click the on the container.

Select **SHOW ROUTING**



MAKE SURE THE **IN-LAB LOCATION** fits with the **DESTINATION** you're trying to set.<sup>17</sup>

<sup>17</sup> For this example, it would not appear if **BH REFERRALS** was set as the destination in **TRANSFER SPECIMENS**.

### I Still Can't Find the Container

If all of the above options fail, try re-logging the sample into the laboratory.

Open Specimen Log-in

Scan the barcode of the sample.

Set the location to your laboratories LOG-IN LOCATION.

Click the Log-in Button.


THE REASON THIS *may* work is because the sample may have been placed on a transfer list that cannot be easily found.


### I Can't Find a List

#### Check the List Options

Many times, this issue is caused because either the **TO LOCATION**, **FROM LOCATION** or **DATE** are incorrectly entered.

Check the **FROM** Location.

 *Check* the **To** Location.

 *Check* the **DATE**.

 *Click* 

Locations					
From:		To:		Transfer Temp:	
From SMC Login		To BH		(All)	
				<input checked="" type="radio"/> All containers <input type="radio"/> Aliquots only <input type="radio"/> Specimens only	
				<input type="button" value="Retrieve"/>	
From Available					
	Accession	Order	Patient Name	Patient ID	
1	15-264-000058A	BUNU	TESTING, NOTWC	1017792	100
2	15-293-000039A	HBA1C	TEST, VON TESTINGSTON	6038342000	3 ml
To Selected					
	Accession	Order	Patient Name		

