


CERNER TRAINING MANUALS

MODIFY COLLECTIONS

MODIFY COLLECTIONS

Modify Collections is the application used to update the collection details of a container.




IT CAN BE USED to update:

- **COLLECTION TIME:**
- **COLLECTION DATE:**
- **COLLECTION ID:**
- **COLLECTION METHOD:**

Using Modify Collections

Enter the ACCESSION NUMBER.

Click 


Double Click on the field to modify.

Enter the new value.

Click  when finished.

IF THE RESULTS have already been verified, Cerner will display a warning.¹

¹This warning is for situations where the REFERENCE RANGE is determined, in part, by the time of collection.



Click .

FINALLY, IT WILL ask if you'd like some labels.

