

CERNER TRAINING MANUALS


TAT MONITOR


CONTENTS


| | | | |
|------------------------|---|----------------------------|---|
| How it Works | 1 | What's Displayed | 2 |
|------------------------|---|----------------------------|---|

TURNAROUND TIME MONITOR

The TAT MONITOR is used to view outstanding orders which have been received in the laboratory.



OPEN TURNAROUND TIME MONITOR by clicking the  icon from the App-Bar.¹

¹  Refer to the APP-BAR PROCEDURE if you need help adding it.




How it Works

Just like PENDING INQUIRY, the TAT MONITOR is pulled using a  SECTION,  SUBSECTION,  BENCH, or  INSTRUMENT. However, in addition, the TAT MONITOR can be pulled using a  SITE.²

² This will include all of the laboratories  SECTIONS and  SUBSECTIONS.

 *Select* a SERVICE RESOURCE in the SERVICE RESOURCE FIELD.



 *Click* 

PathNet Lab Management: TAT Monitor

Task View Help

Selected Service Resource: SSW Lab As Of: 12/3/2015 12:52

| | Received | Name | Accession | Order Mem... | Priority | Location | Dept Status |
|----|------------------|------------------------|---------------|--------------|--------------|----------|-------------|
| 1 | 12/3/2015 11:15 | TEST, VON TESTINGSTON | 15-337-000049 | UAS | RT - Routine | ICUA SMC | In-Lab |
| 2 | 12/3/2015 11:13 | TEST, VON TESTINGSTON | 15-337-000050 | PTT | RT - Routine | ICUA SMC | In-Lab |
| 3 | 10/20/2015 14:20 | TEST, SSW/70707 | 15-293-000062 | ADIFF | ST - Stat | ED SSW | In-Lab |
| 4 | 10/20/2015 14:20 | TEST, SSW/70707 | 15-293-000062 | CBC | ST - Stat | ED SSW | In-Lab |
| 5 | 5/14/2015 11:25 | ZZTESTGENLAB, SMCPC11Y | 15-134-000043 | CBC | ST - Stat | 7N SMC | In-Lab |
| 6 | 5/14/2015 11:25 | ZZTESTGENLAB, SMCPC11Y | 15-134-000043 | ADIFF | ST - Stat | 7N SMC | In-Lab |
| 7 | 12/3/2015 12:25 | Ref. Gen Lab | 15-337-000064 | PHOTIM | ST - Stat | ICUA SMC | In-Lab |
| 8 | 12/3/2015 12:42 | TEST, VON TESTINGSTON | 15-337-000065 | D-Dimer | RT - Routine | ICUA SMC | In-Lab |
| 9 | 12/3/2015 12:51 | TEST, VON TESTINGSTON | 15-337-000064 | ADIFF | RT - Routine | ICUA SMC | In-Lab |
| 10 | 12/3/2015 12:51 | TEST, VON TESTINGSTON | 15-337-000064 | CBC | RT - Routine | ICUA SMC | In-Lab |

Notification Color Legend:

| | Display | Warning | Alert / Expiration |
|----------|---------|---------|--------------------|
| Selected | Display | Warning | Alert / Expiration |
| Normal | Display | Warning | Alert / Expiration |

Ready Refresh: 1 min A3291 TPERR005 12/3/2015 12:53

What's Displayed

The items displayed will be different for each Laboratory. Larger sites will only display STAT and possibly TIME STUDY orders.

Smaller sites will all orders.

